

# Amplite™ Colorimetric Alkaline Phosphatase Assay Kit

## \*Yellow Color\*

### Ordering Information:

Product Number: #11950 (500 assays)

### Storage Conditions:

Keep in freezer and avoid light

### Instrument Platform:

Absorbance microplate readers

## Introduction

Alkaline phosphatase is a highly sensitive enzyme for ELISA, immuno-histochemical, Northern, Southern and Western blot applications. It is widely used in various biological assays (in particular, immunoassays) and ELISA-based diagnostics. This Amplite™ Alkaline Phosphatase Assay Kit uses *p*NPP, a chromogenic phosphatase substrate, to quantify alkaline phosphatase activity in solutions, in cells, as well as on solid surfaces (such as PVDF membranes). The kit provides all the essential components with our optimized “mix and read” assay protocol that is compatible with HTS liquid handling instruments.

This Amplite™ Alkaline Phosphatase Assay Kit can be performed in a convenient 96-well or 384-well microtiter-plate format and easily adapted to automation with no separation steps required. Its signal can be easily read by absorbance microplate reader around 400 nm.

### Kit Key Features

<b>Optimized:</b>	Optimized conditions for detecting alkaline phosphatase activity.
<b>Continuous:</b>	Easily adapted to automation with no separation required.
<b>Convenient:</b>	Formulated to have minimal hands-on time. No wash is required.
<b>Non-Radioactive:</b>	No special requirements for waste treatment.

## Kit Components

Components	Amount
Component A: <i>p</i> NPP (light sensitive)	1 vial
Component B: Assay Buffer	1 bottle (25 mL)
Component C: Alkaline Phosphatase Standard	1 vial (lyophilized powder, 10 units)

## Assay Protocol for One 96-Well Plate

### Brief Summary

**Prepare assay reaction mixture (50 µL) → Add alkaline phosphatase standards or test samples (50 µL) → Incubate at RT or 37°C for 5-30 min → Read absorbance at 400 nm**

*Note: Thaw all the kit components to room temperature before starting your experiment.*

### 1. Prepare *p*NPP stock solutions:

- 1.1 *p*NPP stock solution (100X): Add 300 µL of distilled H<sub>2</sub>O into the vial of *p*NPP (Component A). Mix the reagents well. The stock solution should be used promptly. Any remaining solution need be aliquoted and refrozen at -20°C.

*Note: Avoid repeated freeze-thaw cycles. It will be good for 3-4 weeks if stored at -20°C.*

### 2. Prepare *p*NPP reaction mixture:

- 2.1 Prepare *p*NPP reaction mixture according to the following table and kept from light:

**Table 1.** *p*NPP mixture for one 96-well plate (2X)

Components	Volume
<i>p</i> NPP stock solution (100X, from step 1.1)	50 µL
Assay Buffer (Component B)	5 mL
Total volume	5 mL

*Note: Prepare fresh reaction mixture for each experiment.*

ABD Bioquest, Inc., 923 Thompson Place, Sunnyvale, CA 94085. Tel: 408-733-1055; Fax: 408-733-1304


Ordering: [sales@abdbioquest.com](mailto:sales@abdbioquest.com); 800-990-8053 or 408-733-1055; Fax: 408-733-1304

Technical Support: [support@abdbioquest.com](mailto:support@abdbioquest.com); 408-733-1055



## Data Analysis

The absorbance in blank wells (with equal volume of *p*NPP and H<sub>2</sub>O-0.1%BSA only) is used as a control, and is subtracted from the values for those wells with alkaline phosphatase reactions. The typical data are shown in Figure 1 (alkaline phosphatase standard curve).



**Figure 1.** Alkaline phosphatase dose response on a white/clear bottom 96-well plate using a NovoStar microplate reader (BMG Labtech) measured with the Amplite™ Alkaline Phosphatase Assay Kit. As low as 0.03 mU/well of alkaline phosphatase can be detected with 30 minutes incubation time (n=3).

## References:

1. Zhu X, Jiang C. (2006) 8-Quinoyl phosphate as a substrate for the fluorimetric determination of alkaline phosphatase. *Clin Chim Acta*.
2. Ali AT, Penny CB, Paiker JE, Psaras G, Ikram F, Crowther NJ. (2006) The effect of alkaline phosphatase inhibitors on intracellular lipid accumulation in preadipocytes isolated from human mammary tissue. *Ann Clin Biochem*, 43, 207.
3. Lee DH, Lim BS, Lee YK, Yang HC. (2006) Effects of hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) on alkaline phosphatase activity and matrix mineralization of odontoblast and osteoblast cell lines. *Cell Biol Toxicol*, 22, 39.
4. Ali AT, Penny CB, Paiker JE, van Niekerk C, Smit A, Ferris WF, Crowther NJ. (2005) Alkaline phosphatase is involved in the control of adipogenesis in the murine preadipocyte cell line, 3T3-L1. *Clin Chim Acta*, 354, 101.
5. Rieu JP, Ronzon F, Place C, Dekkiche F, Cross B, Roux B. (2004) Insertion of GPIanchored alkaline phosphatase into supported membranes: a combined AFM and fluorescence microscopy study. *Acta Biochim Pol*, 51, 189.
6. Palermo C, Manduca P, Gazzero E, Foppiani L, Segat D, Barreca A. (2004) Potentiating role of IGFBP-2 on IGF-II-stimulated alkaline phosphatase activity in differentiating osteoblasts. *Am J Physiol Endocrinol Metab*, 286, E648.

**Warning: This kit is only sold for the end users. Neither resale nor transfer to a third party is allowed without written permission from ABD Bioquest. Chemical analysis of kit components is strictly prohibited. Please call us at 408-733-1055 or e-mail us at [info@abdbioquest.com](mailto:info@abdbioquest.com) if you have any questions.**