

# Amplite™ Glucose Quantitation Kit \*Red Fluorescence\*

**Ordering Information:**

Product Number: #40005 (500 assays)

**Instrument Platform:**

Fluorescence microplate readers

**Storage Conditions:**

Keep in freezer and avoid light

## Introduction

Glucose, a monosaccharide, is the most important carbohydrate in biology. It is a source of energy and metabolic intermediate for cell growth. Glucose is one of the main products of photosynthesis and starts cellular respiration in both prokaryotes and eukaryotes. Glucose level is a key diagnostic parameter for many metabolic disorders. This Amplite™ glucose assay kit provides a quick and sensitive method for the measurement of glucose in various biological samples (e.g., serum, plasma, body fluid, food, growth medium, etc.). The assay is robust, and can be readily adapted for high-throughput assays in a wide variety of applications that require the measurement of glucose. For example, the assay might be suitable for monitoring glucose level during fermentation and glucose feeding in protein expression processes. It might also be used for monitoring glucose transporters. In addition, this assay has very low background since it is run in the red visible range that significantly reduces the interference from biological samples. The assay has demonstrated high sensitivity and low interference with 570 nm excitation 590 nm emission. With the Amplite™ Glucose Quantitation Kit, we have detected as little as 3  $\mu$ M D-glucose.

The Amplite™ Glucose Quantitation Kit can be performed in a convenient 96-well or 384-well microtiter-plate format and easily adapted to automation with no separation steps required. The kit uses our Amplite™ HRP substrate that making the kit recordable in a dual mode, the fluorescent signal can be easily read by either fluorescence microplate reader with Ex/Em = 530 to 570 nm/590 to 600 nm (maximum Ex/Em = 540 nm/590 nm) or absorbance microplate reader at 576 $\pm$ 5 nm.

### Kit Key Features

<b>Sensitive:</b>	The kit detects as low as 3 $\mu$ M D-glucose in solution.
<b>Continuous:</b>	Easily adapted to automation with no separation required.
<b>Convenient:</b>	Formulated to have minimal hands-on time. No wash is required.
<b>Non-Radioactive:</b>	No special requirements for waste treatment.

## Kit Components

Component	Amount
Component A: Amplite™ Red (light-sensitive)	1 vial
Component B: Assay Buffer	1 bottle (50 mL)
Component C: Horseradish Peroxidase (HRP)	1 vial (10 units)
Component D: Glucose Oxidase	1 vial (100 units)
Component E: DMSO	1 vial (200 $\mu$ L)
Component F: Glucose	1 vial (144 mg)

ABD Bioquest, Inc., 923 Thompson Place, Sunnyvale, CA 94085. Tel: 408-733-1055; Fax: 408-733-1304

Ordering: [sales@abdbioquest.com](mailto:sales@abdbioquest.com); 800-990-8053 or 408-733-1055

Technical Support: [support@abdbioquest.com](mailto:support@abdbioquest.com); 408-733-1055

## Assay Protocol (for 1 plate)

### Brief Summary

Prepare assay reaction mixture (50  $\mu$ L) → Add Glucose standards or test samples (50  $\mu$ L)  
→ Incubate at 37°C for 10-30 min → Read fluorescence at Ex 540 nm/Em 590 nm

*Note: Thaw all the kit components to room temperature before starting the experiment.*

#### 1. Prepare stock solutions:

- 1.1 **Amplite™ Red stock solution (250X):** Add 100  $\mu$ L of DMSO (Component E) into the vial of Amplite™ Red substrate (component A). The stock solution should be used promptly. Any remaining solution need be aliquoted and refrozen at -20°C.

*Note 1: Avoid repeated freeze-thaw cycles.*

*Note 2: The Amplite™ Red substrate is unstable in the presence of thiols such as dithiothreitol (DTT) and 2-mercaptoethanol. The final concentration of DTT or 2-mercaptoethanol in the reaction should be no higher than 10  $\mu$ M. It is also unstable at high pH (>8.5). The reactions should be performed at pH 7–8. The provided reaction buffer, pH 7.4, is recommended.*

- 1.2 **10 U/ml HRP stock solution:** Add 1 mL of assay buffer (Component B) into the vial of horseradish peroxidase (Component C).

*Note: The unused HRP solution should be divided as single use aliquots and stored at -20°C.*

- 1.3 **100U/ml glucose oxidase solution:** Add 1 mL of assay buffer (Component B) into the vial of glucose oxidase (Component D).

*Note: The unused glucose oxidase solution should be divided as single use aliquots and stored at -20°C*

- 1.4 **400 mM glucose stock solution:** Add 2 mL of assay buffer (Component B) into the vial of glucose (Component F).

*Note: The unused glucose solution can stored at -20°C*

#### 2. Prepare assay reaction mixture:

- 2.1 Prepare Assay reaction mixture according to the following table and kept from light:

**Table 1.** Assay reaction mixture for one 96-well plate (2X)

Components	Volume
Amplite™ Red stock solution (250X, from step 1.1)	20 $\mu$ L
10 U/ml HRP (from step 1.2)	100 $\mu$ L
100 U/ml glucose oxidase (from step 1.3)	100 $\mu$ L
Assay Buffer (Component B)	4.78 mL
Total volume	5 mL

**Table 2.** Layout of glucose standards and test samples in a solid black 96-well microplate:

BL	BL	TS	TS	....	....						
GS1	GS1	....	....	....	....						
GS2	GS2										
GS3	GS3										
GS4	GS4										
GS5	GS5										
GS6	GS6										
GS7	GS7										

*Note: GS= Glucose standards, BL=Blank control, TS=test samples.*

**Table 3.** Reagent composition for each well:

Glucose Standard	Blank Control	Test Sample
Serial dilutions* (50 $\mu$ L)	Assay buffer (Component B): 50 $\mu$ L	50 $\mu$ L

*\*Note 1: Add the serially diluted glucose standards from 3  $\mu$ M to 200  $\mu$ M into wells from GS1 to GS7 in duplicate.*

*Note 2: High concentration of glucose( e. g., 500  $\mu$ M, final concentration) may cause reduced fluorescence signal due to the overoxidation of Amplitude™ red substrate (to a non-fluorescent product).*

### 3. Run Glucose assay

3.1 Prepare a glucose standard by diluting the appropriate amount of the 400 mM glucose stock solution (prepared from step 1.4) into assay buffer (Component B) to produce glucose concentrations of 0 to 200  $\mu$ M, each in a volume of 50  $\mu$ L. A no-glucose buffer control is included as blank control. The final glucose concentrations will be twofold lower (i.e., 0 to 100  $\mu$ M).

3.2 Add 50  $\mu$ L of assay reaction mixture (from step 2) to each well of the glucose standard, blank control, and test samples (see step 2, table 3) so that the total glucose assay volume is 100  $\mu$ L/well

*Note: For a 384-well plate, add 25  $\mu$ L sample, 25  $\mu$ L of assay reaction mixture per well.*

3.3 Incubate the reaction for 10 to 30 minutes at 37°C, protected from light.

3.4 Monitor the fluorescence increase with 530-570 nm (optimal at 540) excitation and 590-600 nm emission using a fluorescence plate reader.

*Note: The contents of the plate can also be transferred to a white clear bottom plate and read by absorbance microplate reader at the wavelength of 576 $\pm$ 5 nm. The absorption detection has lower sensitivity compared to fluorescence reading.*

### 4. Run Data Analysis

The fluorescence in blank wells (with the assay buffer only) is used as a control, and is subtracted from the values for those wells with glucose reactions. The typical data are shown in Figure 1 (glucose standard curve).

*Note: The fluorescence background increases with time, thus it is important to subtract the fluorescence intensity value of the blank wells for each data point.*




Figure 1. Glucose dose response on 96-well black plate using a Novostar microplate reader (BMG Labtech) measured with Amplite™ Glucose Quantitation Kit. As low as 3 µM of glucose can be detected with 30 minutes incubation time (n=3).

## References:

1. Delva P, Degan M, Trettene M, Lechi A. (2006) Insulin and glucose mediate opposite intracellular ionized magnesium variations in human lymphocytes. *J Endocrinol*, 190, 711.
2. Delva P, Degan M, Pastori C, Faccini G, Lechi A. (2002) Glucose-induced alterations of intracellular ionized magnesium in human lymphocytes. *Life Sci*, 71, 2119.
3. Wang XT, Au SW, Lam VM, Engel PC. (2002) Recombinant human glucose-6-phosphate dehydrogenase. Evidence for a rapid-equilibrium random-order mechanism. *Eur J Biochem*, 269, 3417.
4. Leira F, Louzao MC, Vieites JM, Botana LM, Vieytes MR. (2002) Fluorescent microplate cell assay to measure uptake and metabolism of glucose in normal human lung fibroblasts. *Toxicol In Vitro*, 16, 267.
5. Goi G, Bairati C, Burlina A, Massacesi L, Monciotti C, Segalini G, Testa R, Lombardo A. (2000) Plasma glycohydrolase levels in patients with type 1 diabetes at onset and in subjects undergoing an intravenous glucose tolerance test. *Metabolism*, 49, 1352.
6. Manduteanu I, Voinea M, Serban G, Simionescu M. (1999) High glucose induces enhanced monocyte adhesion to valvular endothelial cells via a mechanism involving ICAM-1, VCAM-1 and CD18. *Endothelium*, 6, 315.
7. Ruiz-Munoz LM, Vidal-Vanaclocha F, Lampreabe I. (1997) Enalaprilat inhibits hydrogen peroxide production by murine mesangial cells exposed to high glucose concentrations. *Nephrol Dial Transplant*, 12, 456.